

# HALYARD\*

## Sterilisation Pouches



# MEET YOUR STERILISATION POUCH NEEDS WITH CONFIDENCE

The HALYARD\* portfolio offers a line of high-quality single-use sterilisation pouches and reels to meet your sterile processing department's sterilisation packaging needs.

- Feature triple heat-sealed rail seals
- CE marked and comply with the international norms and standards – ISO 11607 & EN 868
- Allow for aseptic presentation of contents

- Available in:
  - A wide range of sizes
  - Pre-Cut Heat-Seal Pouches
  - Heat Seal Reels (Rolls)
- Feature:
  - **Pre-Cut pouches** feature heat-sealed corners on the thumb notch end, which help prevent dust from collecting at the opening seal and reducing the risk of contamination when pouch is opened
  - **Thumb Notch** for easy opening
  - **Specially designed peel-off seal** that makes opening easy, clean, & fiber-free



# HEAT-SEAL POUCHES AND REELS

HALYARD\* Paper-Film Sterilisation Pouches are made from premium 68 GSM paper consisting entirely of primary virgin fiber with excellent bacterial barrier properties that allow for proper sterilant penetration. The product provides clean peeling performance and high seal strength. Multilayer transparent film provides puncture and abrasion resistance while also delivering a high-quality seal. The combination of our paper and film create a high-performance, durable sterile barrier system that provides users with a reliable clean peel experience. HALYARD\* Paper/Film Pouches are printed with sterilisation process indicators for steam and EO sterilisation modalities.



## HALYARD\* Heat-Seal Paper/Film Pouches

Product Code	Size	Compatible Sterilisation Modality	Packaging
48516	50 x 250mm	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48518	75 x 200mm	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48520	75 x 300mm	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48495	100 x 150mm	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48498	100 x 270mm	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48499	100 x 300mm	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48500	100 x 380mm	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48501	100 x 400mm	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48503	120 x 400mm	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48506	150 x 300mm	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48507	150 x 380mm	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48558	150 x 350mm	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48508	150 x 400mm	STEAM/EO	250 eaches/box 8 boxes/case 2,000 eaches/case
48509	200 x 610mm	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48510	205 x 400mm	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48511	220 x 300mm	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48512	250 x 380mm	STEAM/EO	250 eaches/box 4 boxes/case 1,000 eaches/case
48513	250 x 500mm	STEAM/EO	250 eaches/box 2 boxes/case 500 eaches/case
48514	300 x 500mm	STEAM/EO	250 eaches/box 2 boxes/case 500 eaches/case
48515	420 x 500mm	STEAM/EO	250 eaches/box 2 boxes/case 500 eaches/case

## HALYARD\* Heat-Seal Paper/Film Gusseted Pouches/Reels

48492	150 x 50 x 400mm	STEAM/EO	100 eaches/box 10 boxes/case 1,000 eaches/case
48494	150mm x 50mm x 100m	STEAM/EO	Shrink-wrapped individually 2 reels/case

## HALYARD\* Heat-Seal Paper/Film Reels

Product Code	Size	Compatible Sterilisation Modality	Packaging
48526	75mm x 200m	STEAM/EO	Shrink-wrapped individually 2 reels/case
48521	100mm x 200m	STEAM/EO	Shrink-wrapped individually 1 reels/case
48522	150mm x 200m	STEAM/EO	Shrink-wrapped individually 1 reels/case
48523	200mm x 200m	STEAM/EO	Shrink-wrapped individually 1 reels/case
48524	250mm x 200m	STEAM/EO	Shrink-wrapped individually 1 reels/case
48555	300mm x 200m	STEAM/EO	Shrink-wrapped individually 1 reels/case



Pouch Composition: Paper Backer - Green Tint Film Front