

**Nurse Stethoscope**

REF:	76.001.00.000	black
	76.001.00.001	blue
	76.001.00.002	green
	76.001.00.003	red

Unit: 1 pcs

Packing Unit: 100 pcs

**Description**

Nurse Stethoscope used for artery auscultation during blood pressure measurement

**Material Composition**

Product: stethoscope: Aluminium, stainless steel, PVC (DEHP), Zinc

Packing: individual packing: PE bag, cardboard  
outer case: corrugated cardboard

**Sterilisation**

Non-sterile