



Medical mask, disposable

TD-30-I.2.c-3.2

edition I

date: 2021-05-26

change of day: N/A

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland; T: +48 32 376 07 00-02,06 F: +48 32 370 38 08; www.zarys.com

CLASSIFICATION

Class I medical device compliant with Medical Devices Regulation (EU) 2017/745 (MDR)

INTENDED PURPOSE

- disposable medical device covering the mouth and nose providing a barrier to minimize the direct transmission of infective agents between staff and patient

FEATURES

- made of three layers of non-woven
- modeling nose-piece
- bacteria filtration efficiency BFE \geq 98%
- dimensions: 175 mm x 95 mm
- with earloops or ties
- splash resistant
- packaging: ergonomic cardboard box to remove individual pieces
- available colours: green, blue, white, pink
- type IIR as per: EN 14683
- single-use
- non-sterile
- shelf life: 5 years



CATALOGUE NUMBER	COLOUR	KIND OF MASK	TYPE AS PER EN 14683	SALES PACKAGING BOX	TRANSPORT PACKAGING CARTON
BT-046-R31B	blue	with earloops	type IIR	50 pcs	20 x 50 pcs
BT-046-R31G	green	with earloops	type IIR	50 pcs	20 x 50 pcs
BT-046-R31P	pink	with earloops	type IIR	50 pcs	20 x 50 pcs
BT-046-R31W	white	with earloops	type IIR	50 pcs	20 x 50 pcs
BT-046-R32B	blue	with ties	type IIR	50 pcs	20 x 50 pcs
BT-046-R32G	green	with ties	type IIR	50 pcs	20 x 50 pcs
BT-046-R32P	pink	with ties	type IIR	50 pcs	20 x 50 pcs
BT-046-R32W	white	with ties	type IIR	50 pcs	20 x 50 pcs