

1603833 - Epidurale anesthesieset - 1 st



Specificaties

Merk	Braun
sku	1603833
Inhoud steriele sets	Perican® canule 18G, 1,3 x 80 mm: 1st, Omnifix® spuit 10 ml LS: 1st, Perifix® LS LOR spuit 8 ml: 1st, Omnifix® spuit LS 8 ml: 1st, Omnifix® LS spuit 5 ml: 1st, Omnifix® LS spuit 3 ml: 1st, Sterican® canule 26G, 0,45 x 25 mm, bruin: 1st, Sterican® canule 21G, 0,8 x 40 mm, groen: 1st, Sterican® canule G18, 1,2 x 40 mm, roze: 1st, gaaskompressen 7,5 x 7,5 cm, 8-laags: 5st, swabsticks: 2st, zelfklevende patiëntenlakens, 75 x 90 cm, blauw: 1st, Plateau 20 x 14,8 x 4 cm, groen met web: 1st, wisselplateau rechthoekig, 14,8 x 7,5 cm, 265 ml, wit: 1st, groen veld, 75 x 90 cm: 1st
Europese Regelgeving	MDR - 2017/745/EU - Klasse III
Type verpakking	Set