

**1543939 - PU sloop armkussen - H 15 cm x B 20 cm x L 35 cm - 1
st**



Beschrijving

Afneembare overtrek in katoen en wasbaar in de machine.

Specificaties

| | |
|-------------------------------------|--|
| Merk | Renol |
| SKU leverancier | RCM16099 |
| Normeringen | EN ISO 13485:2016; 14971:2012; 10993:2019; 10993-5:2009; 1041:2008; 15223-1:2016; 62366:2008 |
| Samenstelling | PU |
| sku | 1543939 |
| Kleur sloop | Blauw |
| Materiaal | PU |
| Merk positioneringsmateriaal | Renol |
| Wasbaar | 60° C |
| Europese Regelgeving | MDR - 2017/745/EU - Klasse I |
| Type verpakking | Stuk |