

## 1533616 - Sonicaid FM830 Encore - Acute Fetal/Maternal Monitoring



### Beschrijving

The sonicaid FM830 Encore range is designed to provide flexible solutions for high-risk fetal and maternal assessment.

A high performance fetal and maternal intrapartum monitor uses today's most advanced technology to produce the most accurate acquisition of data for clinician decision support. Whether you need standard fetal monitoring or equipment for specialist, high-risk care, the sonicaid FM830 Encore series gives you accuracy and clarity in every situation.

For close monitoring of the mother as well as the fetus, the Sonicaid FM830 Encore is a fully integrated system which allows the simultaneous monitoring of Maternal ECG, Pulse Oximetry and Blood Pressure without the need for additional standalone devices.

### Specificaties

<b>Merk</b>	Huntleigh
<b>sku</b>	1533616
<b>Type verpakking</b>	Stuk